

Paavan Electronics Industries

PEI Ultrasonic Product Details:



1. Ultrasonic Minicleaner:

- ✓ Power: 100 W
- ✓ Tank capacity- 2.5 L
- ✓ Digital Heater and timer
- ✓ SS Basket with Lid



2. Ultrasonic Bath (6L):

- ✓ Power: 150 W
- ✓ Tank capacity- 6 L
- ✓ Digital Heater and timer
- ✓ SS Basket with Lid
- ✓ Outlet Valve available



3. Ultrasonic Bath (11L):
- ✓ Power: 300 W
 - ✓ Tank capacity- 11 L
 - ✓ Tank Size- 300x250x150mm
 - ✓ Digital Heater and timer
 - ✓ SS Basket with Lid
 - ✓ Outlet Valve available



4. Ultrasonic Bath (24L):
- ✓ Power: 600 W
 - ✓ Tank capacity- 24 L
 - ✓ Digital temp. indicator
 - ✓ Digital Timer
 - ✓ SS Basket with Lid
 - ✓ Outlet Valve available



5. Ultrasonic Bath (30L):
- ✓ Power: 1200 W
 - ✓ Tank capacity- 30 L
 - ✓ Tank Size- 500x300x200mm
 - ✓ Digital temp. Controller
 - ✓ SS Basket with Lid
 - ✓ Outlet Valve available



6. Ultrasonic Bath (60L):

- ✓ Power: 1200 W
- ✓ Tank capacity- 60 L
- ✓ Digital temperature Controller and indicator
- ✓ Outlet Valve available
- ✓ Generator Separately available
- ✓ SS Basket with Lid



Immersible Transducer:

As per the customer requirement, we make immersible transducer by leak proof ultrasonic welding.

Size and material, as per customer requirement.



Ultrasonic Generator:

Ultrasonic MOSFET/ IGBT, compact generator of 1800 W/ 2000 W available with micro-controlled cord for selection of ultrasonic frequency and power display for frequency and timer.