

Since  
20 Years



# PAAVAN ELECTRONICS INDUSTRIES

Ultrasonic Bath  
Sonicator  
Specialist Design  
for Pharma Lab



**We introduced Ourselves a leading manufacturer Of PEI Make Ultrasonic Bath Sonicator of Since 20 Years.**

- Ours is, ISO 9001-2015 company with We have CE Certification For Manufacturing Ultrasonic Bath and Sonicator for Analytical lab and QC lab of Pharma Industries.
- We had Supplied Our PEI Make Ultrasonic Bath Sonicator to Leading pharma Industries in India Of Different Model as Per Their Requirement , We Manufacturer Customised Ultrasonic Bath – Sonicator.



[www.peiultrasonic.co.in](http://www.peiultrasonic.co.in)

## Main Feature of PEI Ultrasonic Bath Sonicator :

- ▶ S.S. Fabrication of Tank and Mfg. Transducers and Ultrasonic Generator in-house facility.
- ▶ We used Very Special Type of Transducers of American Piezo Co. which gives effective Vibration and Work For Longer Duration without Fail .
- ▶ Being Long Experience in Mfg of PEI make Ultrasonic Bath and Technical efficient Service Back Up Our Customer are Always Satisfied .
- ▶ Our Ultrasonic Generator Circuit is Micro Controlled mosfet based with Digital Display for Frequency and Ultrasonic frequency Can Variable Very Feature With latest Technology .

## ULTRASONIC BATH SONICATOR



2.5 Liters, 5 Liters & 10 Liters with Totally S.S. Body. Digital Display for Heater and Timer. With S.S. Basket & LID.



30 Liters Cap with Tank size 500 x 300 x 200 mm and more As per Customer requirement. Totally and S.S. body & Digital Temp Controller and Timer



2.5 Liters, 5 Liters & 10 Liters with Totally S.S. Body. Digital Display for Heater and Timer. With S.S. Basket & LID.



Ultrasonic Bath / Sonicator with cooling By Compressor & Copper condensor Coil Surrounding with SS. Tank, Maintain Temp of Tank up to 5°C, Different size tank available with cooling system as per Customer Requirement Table Top model also available for Sonicator with Cooling



### Ultrasonic Bath Sonicator with Cooling

Ultrasonic Bath / Sonicator with cooling By Compressor & Copper condensor Coil Surrounding with SS. Tank, maintain Temp of Tank up to 5°C, Different size tank available with cooling system as per Customer Requirement Table Top model also available for Sonicator with Cooling



### Ultrasonic Bath Sonicator with Cooling and Shaking

Ultrasonic Bath / Sonicator with Cooling System As earlier describe, Now additional shaking facility also available. Operation of Ultrasonic and Shaking One by one.



**PAAVAN ELECTRONICS  
INDUSTRIES**

**Inquiry : +91 9824132297**

E/54, GIDC Electronic-Estate, Sector-26  
Gandhinagar 382028 (Gujarat)

Email : [peultrasonic@yahoo.co.in](mailto:peultrasonic@yahoo.co.in) • [mdshah2222@gmail.com](mailto:mdshah2222@gmail.com)

**Web : [www.peultrasonic.co.in](http://www.peultrasonic.co.in)**